

REQUIRED INFORMATION

Date of Order: _____ **Quote** ☐ **Order** ☐

Dealer Name: _____

Dealer Account #: _____ PO #: _____

Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

ATP/Therapist: _____

Ship to Address: _____

City: _____ State: _____

Zip: _____


Special Client Conditions: _____

Client Gender: ☐ M ☐ F

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required measurements.

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



Customer Service: 1.800.333.6900 | Fax: 1.800.678.4682 | www.invacare.com

INVACARE POWER BASES

- ☐ **IFX-20R** AVIVA FX Base for Ultra Low Maxx Rehab Seating
HCPCS code K0848 \$7,995
- ☐ **IFX-20SP** AVIVA FX Base for Single Actuator System
HCPCS code K0856 \$7,995
- ☐ **IFX-20MP** AVIVA FX Base for Multiple Actuator System
HCPCS code K0861 \$7,995

Select STF: (Measured from front of seat, with chair level)

☐ **ELOW** 16.75" **LOW** 17.75" **MED** 18.75" **TALL** 19.75" .. N/C

USER WEIGHT LIMITS

- ☐ **U300** Weight Capacity up to 300 lbs STD
- ☐ **UW250** Weight Capacity up to 250 lbs⁽¹⁾ STD

1. Required for Elevate Systems.

PERFORMANCE MOTORS

- ☐ **HSTM** Performance Motors - 75 MPH⁽¹⁾ \$1,400

TRANSPORT TIE DOWN

- ☐ **TRBKTS** Wheelchair Transport Brackets⁽¹⁾ STD
- ☐ **TRRO** Transport Ready Option⁽¹⁾⁽³⁾ \$50
- ☐ **TRRO-E** Transport Ready Option for Elevate⁽¹⁾⁽²⁾⁽³⁾ \$300

1. TRRO Meets ISO-7176 Part 19

2. Required for Elevate Systems to Meet ISO-7176 Part 19, Includes 2 Additional Brackets on Seat Rail.

3. Must include headrest or head array

TIRE OPTION

- ☐ **AB1431-3** 14" x 3" Black Tires w/Gel Foam Inserts N/C

CASTER OPTION

- ☐ **AB1228-3** 8" x 2.5" Black Casters w/Solid Inlay N/C

BATTERY TRAY TYPE

- ☐ **GP24** Group 24 Battery Box N/C

BATTERY CHARGER OPTION

- ☐ **110CHARGER** 110 Volt Battery Charger N/C
- ☐ **L8900** Omnit Battery Charger (\$80)

SHROUD FINISH - MUST PICK ONE

- ☐ **164P**  Standard Black N/C
- ☐ **165P**  Anthracite Grey N/C
- ☐ **166P**  Mercury Blue N/C
- ☐ **167P**  Bengal Red N/C
- ☐ **168P**  Silver Star N/C
- ☐ **169P**  Midnight Purple N/C

RIM INSERT COLORS - MUST PICK ONE

- ☐ **164PR**  Standard Black N/C
- ☐ **165PR**  Anthracite Grey N/C
- ☐ **166PR**  Mercury Blue N/C
- ☐ **167PR**  Bengal Red N/C
- ☐ **168PR**  Silver Star N/C
- ☐ **169PR**  Midnight Purple N/C



ELECTRONICS**NON-EXPANDABLE SINGLE DRIVE JOYSTICKS**

- ☐ **REM110** **LiNX Drive Only Remote**⁽¹⁾⁽²⁾ N/C
(One Drive Function)

- ☐ **NEXC** **Non-expandable Controller** N/C

1. Lights are not available with non expandable remote

2. Standard with rehab seat

NON-EXPANDABLE SINGLE ACTUATOR JOYSTICKS

- ☐ **REM110SS** **LiNX Drive Only Remote**⁽¹⁾ N/C
One drive function with egg switch and mounting hardware to operate single actuator

- ☐ **REM210SS** **LED Non-Expandable Remote**⁽¹⁾ N/C
Three drive functions allow through the drive control operation of single actuator (requires STTDC), includes egg switch and mounting hardware

- ☐ **REM215SS** **LED Non-Expandable Remote with Lights**⁽¹⁾ \$350
Three drive functions allow through the drive control operation of a single actuator (requires STTDC) and lights, includes egg switch and mounting hardware

- ☐ **STTDC** **Function key** HCPCS code E231 \$1,795
for operation of single actuator through the drive control⁽²⁾

1. Lights are not available with non expandable remote

2. Required with REM210SS

NON-EXPANDABLE DUAL ACTUATOR JOYSTICKS

- ☐ **REM110MS** **LiNX Drive Only Remote** N/C
One drive function with 4-way toggle or 4-way push button switch, with mounting hardware, to operate two actuators⁽¹⁾⁽²⁾

- ☐ **REM210MS** **LED Non-Expandable Remote** N/C
Three drive functions allow through the drive control operation of two actuators (requires MTTDC), includes 4-way toggle or 4-way push button switch with mounting hardware⁽²⁾

- ☐ **REM215SS** **LED Non-Expandable Remote with Lights**⁽¹⁾ \$350
Three drive functions allow through the drive control operation of a single actuator (requires STTDC) and lights, includes egg switch and mounting hardware

- ☐ **MTTDC** **Function key** HCPCS code E2311 \$2,995
for operation of two or more actuators through the drive control⁽³⁾

1. Lights are not available with non expandable remote

2. Choose switch or toggle

3. Required with REM210MS

ELECTRONICS - EXPANDABLE

- ☐ **EXPC** **Expandable Controller** HCPCS code E2377 \$700

- ☐ **PWH** **Harness Required for Expandable System**

HCPCS code E2313 \$450

- ☐ **REM211** **LED Expandable Remote**⁽¹⁾⁽²⁾⁽⁶⁾ STD

- ☐ **REM216** **LED Expandable Remote with Lights**⁽¹⁾⁽³⁾⁽⁶⁾ \$350

- ☐ **REM400** **Color 3.5" Touch Screen Remote**⁽⁴⁾⁽⁵⁾⁽⁶⁾ \$1,000

- ☐ **SW400TGL** **Toggle Kit for REM400**⁽⁷⁾ \$200

- ☐ **REM500** **LiNX Color 3.5" Touch Screen**⁽⁴⁾⁽⁵⁾⁽⁶⁾
Display Only (No Driver Control) \$1,500

- ☐ **STTDC** **Function key** HCPCS code E2310 \$1,795
for operation of single actuator through the drive control

- ☐ **MTTDC** **Function key** HCPCS code E2311 \$2,995
for operation of two or more actuators through the drive control

1. Up to 6 actuators, 3 drive functions, not available with specialty drive controls

2. Not available with lights and indicators

3. Must pick lights and indicators

4. Can control lights, must choose

5. 6 Actuators and drive functions

6. Must choose expandable and PWH

7. Adds 2 toggle switches to REM400

ELECTRONICS - EXPANDABLE MODULES

- ☐ **OUT500** **Output Module** \$500

- ☐ **IN200** **Input Module for REM200 Series Remotes** \$500

- ☐ **PKG32666** **Therafin® Sip-N-Puff Breath Tube Kit** \$380

- ☐ **INPUT** **LiNX Input Module & Sip-N-Puff Interface**
(9 Pin Connection) \$1500

- ☐ **COVER-L** **Large Back Shroud Cover**⁽¹⁾ N/C

1. Required when choosing OUT500, IN200 or INPUT options.

LIGHTS AND INDICATORS

- ☐ **LIGHTS** **LED Lights and Indicators**⁽¹⁾ \$700

1. Must order REM216, REM400 or REM500.

JOYSTICK MOUNTING HARDWARE

- ☐ **LEFT** **Left Handed Joystick Mount**⁽¹⁾ N/C

- ☐ **QUAD** **Quad Link Swing-Away height adjustable**

HCPCS code E1028 \$235

- ☐ **RIGHT-Q** **Right Handed Joystick Mount** N/C

- ☐ **LEFT-Q** **Left Handed Joystick Mount** N/C

1. Right handed joystick is standard.

ATTENDANT DRIVE CONTROL CHOICES

- ☐ **ACU** **Proportional Attendant Control** \$950

- ☐ **RIGHT-ACU** **Right Handed Mount** N/C

- ☐ **LEFT-ACU** **Left Handed Mount** N/C

**SPECIALTY PROPORTIONAL CONTROLS
REQUIRES REM400 OR REM500**

- ☐ **CREM** **LiNX Compact Remote** HCPCS code E2373 \$1,500

- ☐ **CREM-LF** **LiNX Compact Remote Low Force**

HCPCS code E2373 \$1,500

**PROPORTIONAL SPECIALTY CONTROL MOUNTING
HARDWARE**

- ☐ **LEFT-AC** **Left Handed Joystick Mount** N/C

- ☐ **RIGHT-AC** **Right Handed Joystick Mount** N/C

LINX ELECTRONIC ACCESSORIES

- ☐ **LAK** **LiNX Access Key**⁽¹⁾ \$100

- ☐ **USBCL** **USB Charger (Mounted in Left Rail)** N/C

- ☐ **USBCR** **USB Charger (Mounted in Right Rail)** N/C

1. Required for programming

JOYSTICK TOPS - HCPCS Code E2323

- ☐ **PC101A** **Bodypoint U Shaped Handle 3"**⁽¹⁾ \$120

- ☐ **PC102A** **Bodypoint U Shaped Handle 4"**⁽¹⁾ \$120

- ☐ **PC107A** **Bodypoint Rubber Dome**⁽¹⁾ \$120

- ☐ **PC110A** **Grooved Mushroom Joystick**

Handle 3/16" Stem⁽¹⁾ \$120

- ☐ **1560** **T Handle Flexible**

Joystick Extension⁽¹⁾ \$120

- ☐ **1561** **Straight Handle Flexible**

Joystick Extension⁽¹⁾ \$120

1. Only available on the REM400 Remote



Please note that the HCPCS codes indicated on this order form are our interpretation of the new codes and how we believe our products should be coded. The individual DMERCs will determine if our interpretation is appropriate and whether or not items coded will be accepted or denied. Motion Concepts cannot guarantee that use of the designated codes will ensure coverage or payment. Contact the insurer to confirm compliance with coding and coverage policy

SYSTEMS MUST SELECT ONE SYSTEM

NOTE: Systems Include Seat Pan, Back Pan, Standard Armrests and Front Rigging. They are designed for use with a 2"-3" cushion.

<input type="checkbox"/>	ULMT	Motion Concepts Ultra Low Maxx Tilt Only (50°) HCPCS code E1002.	\$6,595
<input type="checkbox"/>	ULMR	Ultra Low Maxx Recline Only (168°) with Extended Shear Reduction (ESR) HCPCS code E1004	\$10,795
<input type="checkbox"/>	ULMTR	Ultra Low Maxx Tilt (50°) and Recline (168°) with Extended Shear Reduction (ESR) HCPCS code E1007	\$10,795
<input type="checkbox"/>	AELM	Motion Concepts Add 12" Elevating Lift Seat Module to Tilt and/or Tilt & Recline System⁽¹⁾ HCPCS code E2300	\$2,495

Tilt reduced to 45 degrees with elevate option

1. Up to 250 lbs.

MODULES MUST SELECT SYSTEM-TO ADD MODULES

<input type="checkbox"/>	0SA	0 Degree Tilt Angle	N/C
<input type="checkbox"/>	5FAA	5 Degree Fixed Anterior Assist⁽¹⁾	N/C
<input type="checkbox"/>	5SA	5 Degree Fixed Pre-Tilt⁽²⁾	N/C
<input type="checkbox"/>	8PCM	8 Degree Precline Module⁽³⁾	N/C

1. STF raised 2.2" (with chair tilted to level)

2. STF raised 2" @ front of seat (chair at 5 degree tilt)

3. Reduces maximum recline angle by 8 degrees

ADJUSTABLE SEAT WIDTH RANGE

<input type="checkbox"/>	1620W	16" - 20" Width Range	N/C
<input type="checkbox"/>	2022W	20" - 22" Width Range	N/C

SEAT WIDTH SETTINGS

<input type="checkbox"/>	W16	16" Wide	N/C
<input type="checkbox"/>	W17	17" Wide	N/C
<input type="checkbox"/>	W18	18" Wide	N/C
<input type="checkbox"/>	W19	19" Wide	N/C
<input type="checkbox"/>	W20	20" Wide	N/C
<input type="checkbox"/>	W21	21" Wide⁽¹⁾ HCPCS code K0108	\$329
<input type="checkbox"/>	W22	22" Wide⁽¹⁾ HCPCS code K0108	\$329

1. No charge on K0856, K0861 - only billable on K0848

ADJUSTABLE SEAT DEPTH RANGE

<input type="checkbox"/>	1620D	16" - 20" Seat Depth Range	N/C
<input type="checkbox"/>	2023D	20" - 23" Seat Depth Range	N/C

SEAT DEPTH SETTINGS

<input type="checkbox"/>	D16	16" Deep	N/C
<input type="checkbox"/>	D17	17" Deep	N/C
<input type="checkbox"/>	D18	18" Deep	N/C
<input type="checkbox"/>	D19	19" Deep	N/C
<input type="checkbox"/>	D20	20" Deep	N/C
<input type="checkbox"/>	D21	21" Deep⁽¹⁾ HCPCS code K0108	\$529
<input type="checkbox"/>	D22	22" Deep⁽¹⁾ HCPCS code K0108	\$529
<input type="checkbox"/>	D23	23" Deep⁽¹⁾ HCPCS code K0108	\$529

1. No charge on K0856, K0861 - only billable on K0848

LAP BELTS AND CHEST STRAPS

<input type="checkbox"/>	PBSSB-48	Push Button Style Lap Belt, Length 48" . . STD	
<input type="checkbox"/>	PBSSB-60	Push Button Style Lap Belt, Length 60" . . N/C	
<input type="checkbox"/>	PBSSB-70	Push Button Style Lap Belt, Length 71" . . N/C	
<input type="checkbox"/>	PLB2	Bodypoint® Padded Lap Belt 2 point (67" length)	\$170
<input type="checkbox"/>	PLB4	Bodypoint® Padded Lap Belt 4 point (67" length)	\$215

MATRIX SEAT CUSHION OPTIONS

<input type="checkbox"/>	LC	Matrx Libra 16"w-20"w HCPCS Code E2624 . .	\$485
<input type="checkbox"/>	LC21UC	Matrx Libra 21"w HCPCS Code E2624	\$525
<input type="checkbox"/>	LC22UC	Matrx Libra 22"w HCPCS Code E2625	\$595
<input type="checkbox"/>	PS-V	Matrx PS 16"w - 20" HCPCS Code E2605 . .	\$345
<input type="checkbox"/>	PS-V21UC	Matrx PS 22"w HCPCS Code E2606	\$695
<input type="checkbox"/>	MA-VI	Matrx Vi 16"w - 21"w HCPCS Code E2607 . .	\$425
<input type="checkbox"/>	MA-VI21UC	Matrx Vi 22"w HCPCS Code E2608	\$475
<input type="checkbox"/>	CCOC	On Chair Cushion Credit	-\$90

MATRIX ELAN HEADREST PAD⁽¹⁾⁽²⁾⁽³⁾ HCPCS code E0955

<input type="checkbox"/>	MEHSP-6	Matrx Elan Standard Pad 6"	\$235
<input type="checkbox"/>	MEHSP-10	Matrx Elan Standard Pad 10"	\$235
<input type="checkbox"/>	MEHSP-14	Matrx Elan Standard Pad 14"	\$235
<input type="checkbox"/>	MEHOP-9	Matrx Elan Occipital Pad 9"	\$235
<input type="checkbox"/>	MEHOP-12	Matrx Elan Occipital Pad 12"	\$235
<input type="checkbox"/>	MEH-4P	Matrx Elan 4-Point Pad 11"W x 10"H" . .	\$325
<input type="checkbox"/>	CAH	Concepts Auto Style Headrest	\$350

1. With standard cover

2. Not compatible with sling back (ASBA-2)

3. If mounting on PB Elite must order removable hardware

HEADREST MOUNTING HARDWARE OPTIONS

<input type="checkbox"/>	MEMAH	Elan Multi-Axis Removable HCPCS code E1028	\$220
<input type="checkbox"/>	CMAH	Motion Concepts Multi-Axis Removable HCPCS code E1028	\$195

FINISHED BACK HEIGHTS

Finished Back Heights for Tilt Only Systems⁽¹⁾⁽³⁾ N/C

<input type="checkbox"/>	18"	<input type="checkbox"/>	19"	<input type="checkbox"/>	20"	<input type="checkbox"/>	21"
	TH18		TH19		TH20		TH21
<input type="checkbox"/>	22"	<input type="checkbox"/>	23"	<input type="checkbox"/>	24"	<input type="checkbox"/>	25"
	TH22		TH23		TH24		TH25

Finished Back Heights for Tilt/Recline Systems⁽²⁾⁽³⁾ N/C

<input type="checkbox"/>	20"	<input type="checkbox"/>	21"	<input type="checkbox"/>	22"	<input type="checkbox"/>	23"
	TRH20		TRH21		TRH22		TRH23
<input type="checkbox"/>	24"	<input type="checkbox"/>	25"	<input type="checkbox"/>	26"	<input type="checkbox"/>	27"
	TRH24		TRH25		TRH26		TRH27

1. 4"-5" Gap from bottom of back pan to seat pan.

2. 6"-7" Gap from bottom of back pan to seat pan

3. Finished back height measured from seat pan to top of back rest.

STANDARD TILT BACK CANES

Straight Tilt Only & ASBA-2 Back Canes ⁽¹⁾

Straight Back Cane Height N/C

☐ 18" ☐ 20" ☐ 22" ☐ 24"
BCH18 BCH20 BCH22 BCH24

1. ASBA-2 back cane height = finished back height with sling

Straight Back Cane Angle Setting ⁽¹⁾

Backrest Angle N/C

☐ 75° ☐ 79° ☐ 85° ☐ 90° ☐ 95°
75BA 79BA 85BA 90BA 95BA
☐ 101° ☐ 106° ☐ 112° ☐ 116°
101BA 106BA 112BA 116BA

5° Angled Tilt Only & ASBA-2 Back Canes ⁽¹⁾

Angled Back Cane Height N/C

☐ 20" ☐ 22" ☐ 24"
BCHA20 BCHA22 BCHA24

1. ASBA-2 back cane height = finished back height with sling

Angled Back Cane Angle Setting ⁽¹⁾

Backrest Angle N/C

☐ 81° ☐ 84° ☐ 91° ☐ 96° ☐ 101°
81BA 84BA 91BA 96BA 101BA
☐ 106° ☐ 110° ☐ 117° ☐ 121°
106BA 110BA 117BA 121BA

BACK OPTIONS

☐ RHBK Standard Rehab Back & Cushion ⁽⁴⁾⁽⁵⁾ N/C
☐ SLNGK Sling Back N/C
☐ RPIP Recessed Planar Interface Plate ⁽²⁾⁽³⁾ \$125
☐ OMB Omit Standard Rehab Back & Cushion ⁽¹⁾⁽⁶⁾ -\$75

1. Back post only, only available with tilt system

2. Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.

3. Not available on tilt only or tilt/elevate systems, for use with recline module.

4. ASBA-2 comes with sling, back and seat pan.

5. Not available with ASBA-2

6. Not available when ordering Matrix or high back

UPGRADEABLE BACK OPTIONS

☐ High Back HCPCS Code E2620 \$595

☐ 16" high ☐ 18" high ☐ 20" high
HB16 HB18 HB20

☐ Matrix Elite Back 16-20"w (3" contour) HCPCS Code E2620 . \$625

☐ 14" high ⁽¹⁾ ☐ 16" high ☐ 18" high ☐ 20" high
PBE14 PBE16 PBE18 PBE20

☐ Matrix Elite HD Back ⁽²⁾

21"w HCPCS Code E2620 20"w HCPCS Code E2621 \$795

☐ 16" high ☐ 20" high
PBEH16 PBEH20

☐ Matrix Elite TR 16-20"w HCPCS Code E2620 \$650

☐ 16" high ☐ 18" high ☐ 20" high
PBETR16 PBETR18 PBETR20

UPGRADEABLE BACK OPTIONS (CONTINUED)

☐ Matrix Elite Deep Back 16-20"w (6" contour)

HCPCS Code E2620 \$760

☐ 14" high ☐ 16" high ☐ 18" high
EDB14 EDB16 EDB18

☐ 20" high x 24" wide
EDB20

☐ Matrix Elite Deep HD Back ⁽²⁾

21"w HCPCS Code E2620 20"w HCPCS Code E2621 \$795

☐ 16" high ☐ 18" high ☐ 20" high
EDBHD16 EDBHD18 EDBHD20

☐ Matrix PB (Posture Back) 16-20"w (5" contour)

HCPCS Code E2620 \$625

☐ 16" high ☐ 20" high
MPB-V16 MPB-V20

☐ Matrix PB Deep 16-20"w (7" contour)

HCPCS Code E2615 \$645

☐ 16" high ☐ 20" high
MPB-D16 MPB-D20

☐ BCOC On Chair Matrix Back Credit-\$125

1. Not compatible with 19" or 20" wide seats

2. Only available in 21" or 22" width

TILT ARMRESTS ⁽¹⁾

☐ ADPAH-S Adult Dual Post Adj. 9.75"-12.75"
w/flat pouch STD

☐ ADPAH-M Adult Dual Post Adj. Height 12.75"-15.75"
w/flat pouch STD

☐ ADPAH-T Adult Dual Post Adj. Height 15.75"-18.75"
w/flat pouch STD

☐ ULMFB Ultra Rail Seat Mnt flip Back Cantilever
MAXX Tilt Arm (9-12.5") \$250

☐ HTE-L Horizontal Arm Tube Extension Left. N/C

☐ HTE-R Horizontal Arm Tube Extension Right N/C

1. Height measured from seat pan to top of tube

RECLINE ARMRESTS ⁽¹⁾⁽²⁾

☐ REFB-L Reclining Two Point Adjustable Height
9.5"-13" STD

☐ REFB-T Reclining Two Point Adjustable Height
13"-16" STD

Maxx Style Cane Mounted Height, Width & Angle Adjustable
Cantilever Flip Back Armrest

☐ MXFB-S Straight 9.5"-13" \$250

☐ MXFB-C Curved 12.5"-16" \$250

☐ HTE-L Horizontal Arm Tube Extension Left. N/C

☐ HTE-R Horizontal Arm Tube Extension Right N/C

1. Height measured from seat pan to top of tube.

2. Not available with ASBA-2.

ARM PADS⁽¹⁾

Standard **STD**

Left ☐ **Full (2.25" x 14")** ☐ **Desk (2.25" x 10")**
STAP-RF STAP-RD

Right ☐ **Full (2.25" x 14")** ☐ **Desk (2.25" x 10")**
STAP-LF STAP-LD

Waterfall **\$20**

Left ☐ **Full (14")** WFALL-LF ☐ **Desk (10")** WFALL-LD

Right ☐ **Full (14")** WFALL-RF ☐ **Desk (10")** WFALL-RD

Ergonomic Arm Troughs HCPCS Code E2209 **\$85**

☐ **Left** ERGOL ☐ **Right** ERGOR

Flat Hand Pad⁽²⁾ **\$55**

☐ **Left** FHL-LT ☐ **Right** FHL-RT

Modular Integral Soft Skin Foam **\$32**

Left ☐ **Full (14")** MISK-LF ☐ **Desk (10")** MISK-LD

Right ☐ **Full (14")** MISK-RF ☐ **Desk (10")** MISK-RD

Modular Startex Covered Visco Foam **\$32**

Left ☐ **Full (14")** MSTEX-LF ☐ **Desk (10")** MSTEX-LD

Right ☐ **Full (14")** MSTEX-RF ☐ **Desk (10")** MSTEX-RD

☐ **OAP** **Omit Arm Pads.** **-\$25**

1. Full length pads ship with full length arm tubes (w/extensions);

ergo trough ships with desk length tubes (w/o extension)

2. Must order ergo trough

ARMPAD ACCESSORIES

Elbow Block Left HCPCS Code K0108 **\$195**

☐ **XS(3.25"x4.25")** ☐ **S(3.5"x5.25")**

EBLT-XS

EBLT-S

☐ **M(3.75"x5.75")** ☐ **L(4.25"x6.5")**

EBLT-M

EBLT-L

Elbow Block Hardware Left⁽¹⁾ **N/C**

Mounting hardware extension

☐ **Vertical 3"-Horizontal 5"**

VEXSHEXS-L

☐ **Vertical 3"-Horizontal 7"**

VEXSHEXL-L

☐ **Vertical 5"-Horizontal 5"**

VEXLHEXS-L

☐ **Vertical 5"-Horizontal 7"**

VEXLHEXL-L

Elbow Block Right HCPCS Code K0108 **\$195**

☐ **XS (3.25" x 4.25")** ☐ **S (3.5" x 5.25")**

EBRT-XS

EBRT-S

☐ **M (3.75" x 5.75")** ☐ **L (4.25" x 6.5")**

EBRT-M

EBRT-L

Elbow Block Hardware Right⁽¹⁾ **N/C**

Mounting hardware extension

☐ **Vertical 3"-Horizontal 5"**

VEXSHEXS-R

☐ **Vertical 3"-Horizontal 7"**

VEXSHEXL-R

☐ **Vertical 5"-Horizontal 5"**

VEXLHEXS-R

☐ **Vertical 5"-Horizontal 7"**

VEXLHEXL-R

1. Recommend long horizontal except for seat mounted cantilever arm

ARMPAD ACCESSORIES (CONTINUED)

Multi Axis Upper Extremity Support⁽¹⁾ **\$225**

☐ **Left MASL**

☐ **Right MASR**

Armrest Spacers (pair)⁽²⁾ **\$250**

☐ **1" Offset OAS1**

☐ **1/2" Offset OAS12**

1. Not available with standard pad

2. Only compatible with recline arm, not Maxx Cantilever or dual post tilt arm

MAXX LATERALS

☐ **FLATL** **Left Fixed Mount Lateral Trunk Support⁽¹⁾**

HCPCS code E0956 **\$100**

Select Left Trunk Support Pad Size HCPCS code E0956 **N/C**

☐ **XS (3.25" x 4.25")** L-XS

☐ **S (3.5" x 5.25")** L-S

☐ **M (3.75" x 5.75")** L-M

☐ **L (4.25" x 6.5")** L-L

☐ **FLATR** **Right Fixed Mount Lateral Trunk Support⁽¹⁾**

HCPCS code E0956 **\$100**

Select Right Trunk Support Pad Size **N/C**

☐ **XS (3.25" x 4.25")** R-XS

☐ **S (3.5" x 5.25")** R-S

☐ **M (3.75" x 5.75")** R-M

☐ **L (4.25" x 6.5")** R-L

1. Not compatible with "deep" Matrix backs or sling back

MAXX STYLE SWING AWAY HARDWARE

Swing Away Trunk Support Hardware⁽⁴⁾

☐ **LMSA** **MAXX Style Swing Away Hardware Left⁽¹⁾**

HCPCS code E1028 **\$195**

Select Left Trunk Support Pad Size **\$100**

☐ **XS (3.25" x 4.25")** L-XS-SW

☐ **S (3.5" x 5.25")** L-S-SW

☐ **M (3.75" x 5.75")** L-M-SW

☐ **L (4.25" x 6.5")** L-L-SW

☐ **RMSA** **MAXX Style Swing Away Hardware Right⁽¹⁾**

HCPCS code E0956 **\$195**

Select Right Trunk Support Pad Size **\$100**

☐ **XS (3.25" x 4.25")** R-XS-SW

☐ **S (3.5" x 5.25")** R-S-SW

☐ **M (3.75" x 5.75")** R-M-SW

☐ **L (4.25" x 6.5")** R-L-SW

☐ **ATL2** **Pair - Additional Telescoping Link⁽²⁾⁽³⁾**

(up to 5 1/2" medial offset) **\$150**

1. Up to 2.75" medial offset

2. Recommended for all backs except rehab back, PB and elite dep backs

3. Medial offset depends on back thickness

4. Must choose fixed mount and pads

MATRIX FIXED LATERALS FOR ELITE AND ELITE TR BACKS

☐ **PBLS-L** **Matrix Standard Fixed Left (7.5" x 5")**

HCPCS code E0956 **\$100**

☐ **PBLS-R** **Matrix Standard Fixed Right (7.5" x 5")**

HCPCS code E0956 **\$100**

☐ **PBOL-L** **Matrix Offset Fixed Left** HCPCS code E0956 . . . **\$110**

☐ **PBOL-R** **Matrix Offset Fixed Right** HCPCS code E0956 . **\$110**

MATRIX SWING AWAY LATERALS - FOR ELITE AND ELITE TR BACKS**Swing Away Trunk Support Hardware⁽⁴⁾**

- ☐ SAEHBL **Matrix Style Swing Away Hardware Left⁽¹⁾**
HCPCS code E0128 \$195

Select Left Trunk Support Pad Size \$100

- ☐ XS (3.25" x 4.25") L-XS-SW ☐ S (3.5" x 5.25") L-S-SW
☐ M (3.75" x 5.75") L-M-SW ☐ L (4.25" x 6.5") L-L-SW

- ☐ SAEHBR **Matrix Style Swing Away Hardware Right⁽¹⁾**
HCPCS code E0956 \$195

Select Right Trunk Support Pad Size \$100

- ☐ XS (3.25" x 4.25") R-XS-SW ☐ S (3.5" x 5.25") R-S-SW
☐ M (3.75" x 5.75") R-M-SW ☐ L (4.25" x 6.5") R-L-SW

- ☐ OSAHL **Matrix Offset Fixed Left** HCPCS code E0956 . . \$195
☐ OSAHR **Matrix Offset Fixed Right** HCPCS code E0956 \$195

FIXED HIP SUPPORTS

- ☐ LHSFL **Hip Support Pad with Fixed Hardware Left**
HCPCS code E0956 \$120

Select Left Hip Pad Size HCPCS code E0956 N/C

- ☐ XS (5.5"W x 4"H) LHS-XS ☐ S (7"W x 6"H) LHS-S
☐ M (9"W x 4"H) LHS-M ☐ L (13"W x 4"H) LHS-L

- ☐ LHSFR **Hip Support Pad with Fixed Hardware Right**
HCPCS code E0956 \$120

Select Right Hip Pad Size N/C

- ☐ XS (5.5"W x 4"H) RHS-XS ☐ S (7"W x 6"H) RHS-S
☐ M (9"W x 4"H) RHS-M ☐ L (13"W x 4"H) RHS-L

REMOVABLE HIP SUPPORT HARDWARE**Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports**

- ☐ MQRHL **MAXX Style Quick Release - Left**
HCPCS code E1028 \$215

Select Left Hip Pad Size HCPCS code E0956 \$120

- ☐ XS (5.5"W x 4"H) LHS-XS-SW ☐ S (7"W x 6"H) LHS-S-SW
☐ M (9"W x 4"H) LHS-M-SW ☐ L (13"W x 4"H) LHS-L-SW

- ☐ MQRHR **MAXX Style Quick Release - Right**
HCPCS code E1028 \$215

Select Right Hip Pad Size HCPCS code E0956 \$120

- ☐ XS (5.5"W x 4"H) RHS-XS-SW ☐ S (7"W x 6"H) RHS-S-SW
☐ M (9"W x 4"H) RHS-M-SW ☐ L (13"W x 4"H) RHS-L-SW


STANDARD FRONT RIGGING

- ☐ **Basic Fixed Center Mount Foot Platform⁽¹⁾⁽²⁾** N/C

☐ 70° FCMT70

☐ 90° FCMT90

☐ 97° FCMT97



- ☐ FCP **Flat Calf Pads for Fixed Center Mount Foot Platform** \$350

1. Comes standard with 12" x 11" footplate (no heel loops)

2. Seat pan to footplate 9" - 17"

POWER CENTER MOUNT FOOT PLATFORM OPTIONS

- ☐ LNX **Power Elevating Center Mount Foot Platform⁽¹⁾⁽²⁾** HCPCS code E1012 \$2,995

- ☐ **Individual Foot Plate Options for LNX** \$330

☐ Small (5"x7.5") IFPS ☐ Medium (5.5"x9.5") IFPM

☐ Large (6"x11.5")⁽³⁾ IFPL


- ☐ **Heel Loop - Pair** HCPCS code E0951

- ☐ **Foot Platform Options for LNX**

☐ Small (12"W x 8.5" D) LNXSFP \$180

☐ Medium (12"W x 10.75" D) LNXMFP \$180

☐ Extra Large (17"W x 13" D)⁽⁴⁾ LNXELFP \$250



1. KTH range 13"-19"; comes standard with 12"x9" foot platform (no heel loops) and independent calf pads.

2. Not available with 16" seat depth

3. Foot plates may interfere with casters in certain configurations

4. May need to adjust leg position to avoid interference with front casters

MANUAL CENTER MOUNT FRONT RIGGING

- ☐ **IVC Center Mount Adjustable Knee Angle Foot Platform⁽¹⁾**
HCPCS code E0990 \$595

Select Seat Pan to Footboard Length N/C

☐ 9.25" - 12" CTMT-FJ ☐ (12.25" - 15") CTMT-F

Select Individual Foot Plates or Footboard N/C

☐ Small (5.25"x8")

SML

☐ Medium (5.25"x10")

MDM

☐ Large (5.9"x12")

LRG

☐ **Center Mount Footboard**

CTMTFB

1. Not available with 16" seat depth

POWER FRONT RIGGING OPTIONS

- ☐ **MPELR** Maxx Swing-Away Power Elevating
Legrests HCPCS code E1010 \$2,490

Select seat pan to footplate length:

- ☐ S (10"-14") ☐ M (13"-17") ☐ L (16"-20")
MSKTHL-S MSKTHL-M MSKTHL-L

MANUAL FRONT RIGGING OPTIONS

- ☐ **MELR** Maxx Swing-Away Manual Elevating
- Straight HCPCS code E0990 \$680

Select seat pan to footplate length:

- ☐ S (10"-14") ☐ M (13"-17") ☐ L (16"-20")
MSKTHL-S MSKTHL-M MSKTHL-L

SWING AWAY FRONT RIGGING OPTIONS

- ☐ **M70HD** 70° HD Swing Away Footrests⁽¹⁾
(Length 13"-17") \$495

- ☐ **70NHD** 70° Swing Away Footrests⁽¹⁾ N/C

1. Must choose footplate

FOOT PLATE OPTIONS FOR ELEVATING & SWING AWAY FRONT RIGGING

- ☐ **Adjustable Angle Foot Plates (sold as pairs)**
HCPCS code K0040. \$330

- ☐ M (5.25"W x 8"D) ☐ L (6.25"W x 8"D)
MAP LAP

- ☐ **Multi-Axis Adjustable Angle Foot Plates (sold as pairs)**
HCPCS code K0040. \$450

- ☐ M (5.25"W x 8"D) ☐ L (6.25"W x 8"D)
MMA LMA

- ☐ **PHL** Heel Loop HCPCS code K0951. \$40

CALF PANEL & STRAP OPTIONS

- ☐ **CPMD** Medium Bodypoint® Aeromesh™ Calf Strap
(35.5"x9") for chairs 15"-18"W⁽¹⁾ \$115

- ☐ **CPLG** Large Bodypoint® Aeromesh™ Calf Strap
(38.5"x9") for chairs 18"-22"W⁽¹⁾ \$115

- ☐ **PHL** Heel Loop HCPCS code K0951. \$40

1. Only compatible with HD Swing Away Footrest

POWERED SEATING SWITCHES - SINGLE ACTUATOR

- ☐ **2PB** Dual Push Button \$150
☐ **PBMT** Push Button Mounting Kit N/C
☐ **ASL300A** Egg Switch. \$156
☐ **2SW** Dual Toggle. \$250
☐ **ASL611** Arm Rest Switch Mounting Bracket \$252

POWERED SEATING SWITCHES - MULTIPLE ACTUATOR

- ☐ **FWTM** 4-Way Toggle Switch N/C
☐ **QPBM** Quad Push Button Toggle Switch N/C
☐ **10WAY** 10-Way Switch (5 actuator functions)⁽¹⁾ . . \$255
☐ **LEFT-S** Left Handed Mount N/C
☐ **RIGHT-S** Right Handed Mount N/C
☐ **AMSA** Attendant mount on back shroud⁽²⁾ N/C

1. Requires expandable electronics

2. 10-Way switch only

ELECTRONICS OPTIONS

- ☐ **ASL300A** Egg Switch. \$156
☐ **ASL304A** Wobble Switch \$180
☐ **ASL314** Ultralight Mechanical \$180
☐ **2PB** Dual Push Button \$150
☐ **2SW** Dual Toggle. \$250
☐ **MPY** Stereo/Mono Splitter \$15

ADDITIONAL SWITCH MOUNTING HARDWARE

- ☐ **PBMT** Push Button Mounting Kit N/C
☐ **ASL611** Arm Rest Switch Mounting Bracket . . . \$252

VENT TRAYS⁽¹⁾⁽²⁾⁽³⁾

- ☐ **VENT** Articulating Compact Vent Tray \$1,350

1. No 3rd battery or stp down converter available

2. Available on tilt and recline only.

OFF CHAIR OPTIONS

- ☐ **PASL930L** Swing Away Head Rest Mount
w/Egg Switch - Left HCPCS code E1028 \$486
☐ **PASL930R** Swing Away Head Rest Mount
w/Egg Switch - Right HCPCS code E1028 . . . \$486
☐ **PASL930B** Left & Right Swing Away Adjustable Switch
Mount on Head Rest HCPCS code E1028 . . . \$720

MILITARY EMBROIDED PATCHED

- ☐ **TAGPA6** Patch Embroided Military Army \$19.86
☐ **TAGPN6** Patch Embroided Military Navy \$19.86
☐ **TAGPAF6** Patch Embroided Military Air Force . . . \$19.86
☐ **TAGPM6** Patch Embroided Military Marines . . . \$19.86

Total Retail Value \$ _____

NOTES:

[illegible]

NOTES:

[illegible]